

Hotpack Undercounter, Freestanding, and SpaceSaver™ Glassware Washer Accessories

Handling and washing almost any type of laboratory glassware is fast and simple with this broad range of baskets and inserts. Most baskets are constructed of rugged, corrosion-resistant stainless steel. Conversion to spindle washing in spindle ready models is quick and simple using the spindle accessories.



Spindle Rack

Part Number	Description
UCW-155	Lower spindle rack, 22 spindles, 10-1/4" (260 mm) clearance, must be used with UCW-180 or 180H
UCW-180	Upper spindle rack, 22 spindles, 7-3/4" (196 mm) clearance, must be used with UCW-155
UCW-181	Lower spindle rack, 22 spindles, 20-3/4" (527 mm) clearance, must be used without an upper rack
UCW-180H	Half depth upper spindle rack, 12 spindles, 7-3/4" (196 mm) clearance, must be used with UCW-155. Allows washing of short and tall items in the same load
UCW-220	Spindle rack set, includes UCW-155, UCW-180, and parts to convert UCW-180 into UCW-181



Open Basket with Handles and Retainer Top

Part Number	Description
0601-864	Open Basket w/handles, stainless steel, fits 1/2 of top or bottom rack



Petri Dish Insert

Part Number	Description
0601-810	Petri Dish Insert, stainless steel, fits 1/2 of top or bottom rack




Test Tube Baskets and Lids

Part Number	Description
0601-857	Test Tube Basket & lid, up to 150 tubes, 75 mm long, stainless steel, fits 1/4 of top or bottom rack
0601-858	Test Tube Basket & lid, up to 150 tubes, 105 mm long, stainless steel, fits 1/4 of top or bottom rack
0601-859	Test Tube Basket & lid, up to 150 tubes, 165 mm long, stainless steel, fits 1/4 of bottom rack only
0601-860	Test Tube Basket & lid, up to 150 tubes, 200 mm long, stainless steel, fits 1/4 of bottom rack only



Arrow Clip Insert

Part Number	Description
0601-850	Arrow Clip Insert, 14 small & 14 large clips, stainless steel, fits 1/2 of bottom rack only




Small Basket and Cover

Part Number	Description
0601-712	Small Basket & Cover, 7.5" (190 mm) L x 7.5" (190 mm) W x 4" (100 mm) D, stainless steel, fits 1/4 of top or bottom rack



Bottle Insert

Part Number	Description
0601-879	Bottle Insert, full rack, holds 16 1-liter bottles, stainless steel, for use in bottom rack only
0601-880	Bottle Insert, full rack, holds 9 2-liter bottles, stainless steel, for use in bottom rack only
0601-881	Bottle Insert, full rack, holds 48 50 mL bottles, stainless steel, for use in bottom rack only
0601-882	Bottle Insert, full rack, holds 44 100 mL bottles, stainless steel, for use in top or bottom rack
0601-883	Bottle Insert, full rack, holds 24 250 mL bottles, stainless steel, for use in top or bottom rack
0601-884	Bottle Insert, full rack, holds 20 500 mL bottles, stainless steel, for use in top or bottom rack



Other Available Items

Part Number	Description
0601-805	Erlenmeyer Flask Insert, holds 5 flasks up to 1500 mL, use on top or bottom
0601-874	Watch glass Insert, up to 40 pieces, up to 125 mm dia., fits 1/2 of top or bottom rack
0601-807	BOD Bottle Insert, holds 18 bottles, fits 1/2 of top or bottom rack
0601-886	Cover/mesh bottom for full size inserts, stainless steel, for use in top or bottom rack
0601-887	Bottom mesh for 1/2 size inserts, stainless steel, for use in top or bottom rack



Insert for Small Graduate Cylinders

Part Number	Description
0601-815	Graduate Cylinder Insert, holds 32 small cylinders, stainless steel, fits 1/2 of top or bottom rack
0601-904	Retainer top for 0601-807, 810, 815, 864 inserts, stainless steel, for use in top or bottom racks

Field Installed Accessories

Part Number	Description
UCW-DIB	Booster pump for low pressure DI water supply (unit mounted, field installed)
UCW-DNS	Automatic dispensing system for liquid detergent and neutralizer, adjustable with audible alarm for low additive level, 115V



SP Industries Corporate Headquarters
 935 Mearns Road • Warminster, PA 18974 USA
 Toll Free: 800.523.2327 • Phone: 215.672.7800 • Fax: 215.672.7807
 e-mail: cs@SPindustries.com • www.SPindustries.com

VirTis

FTS

Wilmad | LabGlass

Genevac

Hull