



*VirTis Bench Top Freeze Dryers*

*Helping you make more time for science*

**FREE OFFER**

**\$1000**

**in FREE Flaskware**

**Purchase a VirTis Bench Top Manifold Freeze Dryer with your choice of manifold and vacuum pump between now and March 31, 2010 and receive **\$1000 of free flaskware** as our *Thank You*.**



***VirTis Bench Top Freeze Dryers Provide***

***Superior corrosion resistance and vacuum pump protection with condenser temperatures as low as -105C***

***Acrylic manifolds for processing inorganic solvents  
Stainless steel manifolds for processing organic solvents***

***Avoid costly vacuum pump over-sizing and maintenance.***

VirTis freeze dryers and valves are designed with as few vacuum connections as possible to minimize potential leak spots. Avoid vapors by-passing the chilled condenser and prematurely contaminating vacuum pump oil by utilizing a "right sized" vacuum pump,

***Faster turn-around time with hot gas defrost.***

Quick and easy hot gas defrost, less fuss & muss and less waiting for that large block of ice in the condenser to melt.

Flasks available in 40 ml, 75 ml, 150 ml, 300 ml, 600 ml, 1200 ml and 2000 ml borosilicate with optional safety coating.

Interested in ampoules or other freeze drying glassware? Check with your sales representative for a complete listing of freeze drying glassware offered.

**Redeemable at time of system purchase.**



For more information please contact your local VirTis representative  
or  
contact Customer Service at [CSGardiner@SPIndustries.com](mailto:CSGardiner@SPIndustries.com)  
845-255-5000  
800-431-8232